

Syringes Product Catalogue

SHANGHAI MEKON MEDICAL DEVICES CO., LTD

Tel:+008613918384200

Email:Ella@mekonneedle.com

Add:526No 697-3 Lingshi road Shanghai 200072 China

Web:www.mekonmedical.com

SHANGHAI MEKON MEDICAL DEVICES CO., LTD

### **Sterile Syringes for Single Use (3-Parts)**

#### **Structure and Composition:**

The Syringes are assembled by Barrel, Plunger, Piston with/without Hypodermic Needles. All parts and material for this product meet medical requirements. Sterilized by EO.

#### Main Material:

PP, Isoprene Rubber, Silicone Oil

#### Intend Use:

Sterile syringes are intended to inject drug for patients. And syringes are intended for use immediately after filling and are not intended to contain the medicament for extended periods of time.

#### **Certificate and Quality Assurance:**

CE, FDA, ISO13485



Sterile Syringes for Single Use (Luer Slip)

#### **Specification:**

Luer Slip, Luer Lock 1ml, 2ml, 2.5ml, 3ml, 5ml, 10ml, 20ml, 25ml, 30ml, 35ml, 50ml, 60ml, 100ml

#### Refer Code:

Specification	Luer Lock	Luer Slip	Eccentric Tip
1ml	SR11101	SR11201	
2ml	SR11102	SR11202	
2.5ml	SR1112.5	SR1122.5	
3ml	SR11103	SR11203	
5ml	SR11105	SR11205	
10ml	SR11110	SR11210	SR11310
20ml	SR11120	SR11220	SR11320
30ml	SR11130	SR11230	SR11330
35ml	SR11135	SR11235	SR11335
50ml	SR11150	SR11250	SR11350
60ml	SR11160	SR11260	SR11360
100ml	SR111100	SR112100	SR113100



Suitable to be connected with power-driven syringe pumps



Sterile Syringes for Single Use (Luer Lock)

01

### **Sterile Syringes for Single Use (3-Parts)**

#### **Structure and Composition:**

The Syringes are assembled by Barrel, Plunger, Piston with/without Hypodermic Needles. All parts and material for this product meet medical requirements. Sterilized by EO.

#### Main Material:

PP, Isoprene Rubber, Silicone Oil, SUS 304 Stainless Steel Cannula

#### **Intend Use:**

Sterile syringes are intended to inject drug for patients. And syringes are intended for use immediately after filling and are not intended to contain the medicament for extended periods of time.

## **Certificate and Quality Assurance:**

CE, FDA, ISO13485



#### **Specification:**

Luer Slip, Luer Lock 1ml, 2ml, 2.5ml, 3ml, 5ml, 10ml, 20ml, 25ml, 30ml, 35ml, 50ml, 60ml, 100ml

#### **Needle Size:**

15G-33G

Sterile Syringes for Single Use (with Hypodermic Needle)

## Sterile Syringes for Single Use (Fixed Needles)


#### **Structure and Composition:**

The Syringes are assembled by Barrel, Plunger, Piston, Fixed Needles.

All parts and material for this product meet medical requirements. Sterilized by EO.

#### **Certificate and Quality Assurance:**

CE, FDA, ISO13485



## Sterile Safe Syringes for Single Use (Shielded-Needle)


#### **Structure and Composition:**

The Syringes are assembled by Barrel, Plunger, Piston, Safety Needles.

All parts and material for this product meet medical requirements. Sterilized by EO.

#### **Main Material:**

PP, Isoprene Rubber, Silicone Oil, SUS304 Stainless Steel Cannula



#### **Intend Use:**

Sterile syringes are intended to inject drug for patients. And syringes are intended for use immediately after filling and are not intended to contain the medicament for extended periods of time. The attached needle safety shield can be manually activated to cover the needle immediately after use to minimize risk of accidental needle sticks.

## Sterile Safe Syringes for Single Use (Shielded-Needle)

#### **Structure and Composition:**

The Syringes are assembled by Sliding Sleeve, Needle Protective Cap, Needle Tube, Barrel, Plunge and Piston. All parts and material for this product meet medical requirements. Sterilized by EO.



are not intended to contain the medicament for extended periods of time.

The attached needle safety shield can be manually activated to cover the

needle immediately after use to minimize risk of accidental needle sticks.

### **Syringes for Pre-Filled Use**

#### **Structure and Composition:**

The Syringes are assembled by barrel, plunger, piston and syringe cap. All parts and material for this product meet medical requirements.

#### **Main Material:**

PP, BIIR, Silicone Oil

(We can also develop this type of product based on other customized material.)



## Sterile Syringes for Single Use (Colorful Plungers)

#### **Structure and Composition:**

The Syringes are assembled by Barrel, Plunge, Piston with/without Hypodermic Needles. All parts and material for this product meet medical requirements. Sterilized by EO.

#### Main Material:

PP, Isoprene Rubber, Silicone Oil

#### **Specification:**

Luer Slip, Luer Lock 1ml, 2ml, 3ml, 5ml, 10ml, 20ml, 30ml, 50ml, 60ml

## **Certificate and Quality Assurance:**

CE, ISO13485

#### Refer Code:

Specification	Luer Lock	Luer Slip	Eccentric Tip
1ml	SR11101	SR11201	
2ml	SR11102	SR11202	
2.5ml	SR1112.5	SR1122.5	
3ml	SR11103	SR11203	
5ml	SR11105	SR11205	
10ml	SR11110	SR11210	SR11310
20ml	SR11120	SR11220	SR11320
30ml	SR11130	SR11230	SR11330
35ml	SR11135	SR11235	SR11335
50ml	SR11150	SR11250	SR11350
60ml	SR11160	SR11260	SR11360
100ml	SR111100	SR112100	SR113100



#### **Intend Use:**

Sterile syringes are intended to inject drug for patients. And syringes are intended for use immediately after filling and are not intended to contain the medicament for extended periods of time.

## Sterile Syringes for Insulin for Single Use

#### **Structure and Composition:**

The Syringes are assembled by Needle Protective Cap, Needle Tube, Barrel, Plunge, Piston and Protective Cap. All parts and material for this product meet medical requirements. Sterilized by EO.

#### Main Material:

PP, Isoprene Rubber, Silicone Oil, SUS304 Stainless Steel Cannula

#### **Certificate and Quality Assurance:**

Specification: 1ml, 0.5ml, 0.3ml U-40, U-100

CE, FDA, ISO13485

Needle Size:

27G-31G



#### **Intend Use:**

The sterile syringes for insulin for single use intended to be used to administer an injection of insulin to a patient subcutaneously.

#### Refer Code:

Specification	U-40	U-100
0.3ml	INS111	INS112
0.5ml	INS121	INS122
1ml	INS131	INS132

## Sterile Syringes for Insulin for Single Use (Retractable)

#### **Structure and Composition:**

The Syringes are assembled by Sliding Sleeve, Needle Protective Cap, Needle Tube, Barrel, Plunge and Piston. All parts and material for this product meet medical requirements. Sterilized by EO.

#### Main Material:

PP, Isoprene Rubber, Silicone Oil, SUS304 Stainless Steel Cannula

#### **Specification:**

1ml, 0.5ml, 0.3ml U-40, U-100

#### Refer Code:

U-40	U-100	
INS211	1	
INS221	INS222	
INS231	INS232	
	INS211 INS221	

#### **Needle Size:**

27G-31G



#### **Certificate and Quality Assurance:**

CE, FDA, ISO13485

#### Intend Use:

The sterile syringes for insulin for single use intended to be used to administer an injection of insulin to a patient subcutaneously.

09

## Sterile Syringes for Single Use (PC)

#### **Structure and Composition:**

The Syringes are assembled by Barrel, Plunge, Piston with/without Hypodermic Needles. All parts and material for this product meet medical requirements. Sterilized by EO.

#### Main Material:

PC, Isoprene Rubber, Silicone Oil

#### **Specification:**

Luer Lock

1ml, 3ml, 5ml, 10ml, 20ml, 30ml

#### Refer Code:

Specification	Luer Lock
1ml	SR12101
3ml	SR12103
5ml	SR12105
10ml	SR12110
20ml	SR12120
30ml	SR12130



### CE, ISO13485, FDA

Intend Use:

Sterile syringes with needle are intended to inject drug for patient. And syringes are intended to use immediately after filling and are not intended to contain the medicament for extended periods of time.

## Sterile Irrigation Syringes for Single Use

#### **Structure and Composition:**

The Irrigation Syringes are assembled by Barrel, Piston and Plunger. All parts and material for this product meet medical requirements. Sterilized by EO.

#### Main Material:

PP, Isoprene Rubber, Silicone Oil, TPE

#### **Specification:**

Nominal capacity of A type (Pull ring type): 30ml, 50ml, 60ml, 100ml Nominal capacity of B type (Push type): 30ml, 50ml, 60ml, 100ml Nominal capacity of C type (Ball capsule type): 30ml, 50ml, 60ml, 100ml

#### Refer Code:

A type	ISR1
B type	ISR2
C type	ISR3



#### **Certificate and Quality Assurance:**

CE, ISO13485, FDA

#### Intend Use:

Irrigation syringes is used to irrigate a body orifice or wound for transient use, the solution after irrigation will flow out instead of entering the human body.

11 KUL www.kdlchina.com www.kdlchina.com KUL 12

## Sterile Syringes for Single Use (2-Parts)

#### **Structure and Composition:**

The Syringes are assembled by Barrel, Plunge, with/without Hypodermic Needles. All parts and material for this product meet medical requirements. Sterilized by EO.

#### **Main Material:**

PP, Silicone Oil, SUS304 Stainless Steel Cannula

#### **Specification:**

Luer Slip 1ml, 2ml, 3ml, 5ml, 10ml, 20ml

#### **Refer Code:**

Specification	Luer Lock	
1ml	SR21201	
2ml	SR21202	
3ml	SR21203	
5ml	SR21205	
10ml	SR21210	
20ml	SR21220	



#### **Intend Use:**

Sterile syringes with needle are intended to inject drug for patient. And syringes are intended to use immediately after filling and are not intended to contain the medicament for extended periods of time.

### **Sterile Syringes for Veterinary Use**

#### **Structure and Composition:**

The Syringes are assembled by Barrel, Plunger, Piston, Protective Cap, with/without Hypodermic Needles. All parts and material for this product meet medical requirements. Sterilized by EO.

#### **Main Material:**

PP, Isoprene Rubber, Silicone Oil, sus304 Stainless Steel Cannula

#### **Specification:**

With Needles: 3ml, 5ml, 10ml Without Needels: 3ml, 5ml, 10ml, 20ml, 30ml, 60ml

# Specification Luer Lock 3ml VSR103 5ml VSR105 10ml VSR110 20ml VSR120 30ml VSR130 60ml VSR160

Refer Code:



## **Veterinary Needle**

## Specification: Plastic hub

30G	
27G	
26G	
25G	Effective Length
24G	Effective Length 13≤L≤50mm
23G	13363011111
22G	
21G	
20G	

#### Aluminum Hub

18G	Effective Length
16G	Effective Length 13≤L≤50mm
14G	132623011111



13 KDL\* www.kdlchina.com

#### **Aesthetic Products**

#### **Sterile Aesthetic Cannula**

Blunt end, side hole, with scale

#### **Sterile Syringe**

Luer lock, slender diameter, high injection pressure, suitable for injection of high viscosity filler. The finger is designed to fit the large hand peg, and the annular push rod is easy to operate with one hand, which can effectively control the injection speed. High permeability medical materials with good drug compatibility.

#### **Hypodermic Needle**

#### **Certificate and Quality Assurance:**

Microneedle, Size: 30G-34G

CE, ISO13485, FDA

#### Main material

PP / PC, Isoprene Rubber, Silicone oil, SUS304 Stainless Steel Cannula

#### Intend use

These products are intended for hypodermic injection of body fluids or tissues, drugs or medical devices approved for human use.



### Three Ring Sterile syringes

#### Main material

PP, Isoprene Rubber, Silicone oil, ABS

#### **Certificate and Quality Assurance:**

CE, ISO13485

#### Intend use

Sterile syringes are intended to inject drug for patients.

And syringes are intended for use immediately after filling and are not intended to contain the medicament for extended periods of time.

Luer lock, suitable for injection of high viscosity filler.

The finger is designed to fit the large hand peg, and the annula push rod is easy to operate with one hand, which can effectively control the injection speed.



