



CHG APPLICATOR

Product Description

CHG Applicators have been applied for skin antiseptics and Surgery, they are made of medical grade 100PPI open-cell polyurethane foam & ABS handle. It is with formulation of 2% chlorhexidine gluconate (CHG) and 70% isopropyl alcohol (IPA) to meet the strict clinical requirements of FDA, thus it would kill more bacteria than traditional Iodophor & Alcohol Applicators. They have been widely used for Hospitals, INS (Infusion Nurses Society), SHEA (Society for Healthcare Epidemiology of America) and IDSA (Infectious Diseases Society of America). Customization of solution formulation is welcome.

Product Features

- **Rapid-acting:** Chlorhexidine holds better immediate antimicrobial activity.
- **Persistent:** CHG Applicators keep antimicrobial activity for at least 48 hours.
- **Broad spectrum:** MRSA, VRE, CRE, Clostridium Difficile and Acinetobacter.
- **Aseptic technique:** The Applicators reduce the risk of direct hand-to-patient contact.
- **Active in protein-rich biomaterials:** Remains active in the presence of blood and serum.
- **Friction scrub:** Applicators promote gentle friction scrub to make CHG solution penetrate the first five layers of the stratum corneum where 80% of skin-dwelling microorganisms reside.





CHG APPLICATOR

Product Specifications

Item No.	Volume	Total Length	Head Length	Head Width	Head Thickness	Diameter 1	Diameter 2
MCA-105	10.5 ml	185 mm	45 mm	45 mm	8 mm	23 mm	17.5 mm
Solution formulation, 2% Chlorhexidine gluconate and 70% Isopropyl alcohol							
Packing information, 40pcs / case, individual blister packaging, the case dimension is 510 ×490 ×280mm							

Sterilization Process

We take Gamma Ray sterilization process for CHG Applicators if customers have specific requirement.



Miraclean Technology Co.,Ltd.

Address: Room 301, Building A, No.18, Rongshuxia Industrial Zone, Tongxin Community, Baolong Street, Longgang District, Shenzhen, China, 518116

Tel: +86-0755- 89616773 Fax: +86-755-89616775

Web: www.mantacc.com



Share Info Consultant Service LLC Repräsentanzbüro

Address: Heerdter Lohweg 83, 40549 Düsseldorf

Tel: 0049 1767 0057022

Fax: 0049 1767 0057022

Emai: info@share-info.cn

DIMDI NO.: DE/0000047946

CE Logo

