




# Disposable Sampling Swab

## Product Description

93050H is a Foam Sampling Swab used for cell & virus specimen collection, and it has been designed for oropharyngeal clinical diagnostics. This specialty applicator is tipped with 100ppi medical grade polyurethane foam. It is bonded securely to a hot-low PP handle that snaps when bent at a 90° angle. It has been well recognized and adopted by the diagnostic test kit manufacturers who produce reagents in molecular genetics, forensics and clinical laboratories sectors.

## Product Features

- Medical grade quality
- Excellent collection and release properties
- Soft foam tip for patient comfort
- Used dry or with media
- Made in China



## Product Specifications

Item No.	Dimension of foam tip			Dimension of PP handle			
	Width	Thickness	Length	Diameter 1	Diameter 2	Breakpoint	Total length
93050H	6±0.5mm	5.5±0.5mm	17±3mm	3 mm	Taper to 3.0	53/86±3mm	163±15mm
Packaging				100pcs/bag, 50bags/ctn, 5000pcs/ctn			
Sterilization method				Ethylene oxide sterilization			
Ethylene oxide residue				≤10µg/g			
Best storage conditions				4℃— 30℃			
Quality guarantee period				3 years			



# Disposable Sampling Swab

---

## Directions

Step 1: Put the tip of swab into mouth.

Step 2: Squeeze with your tongue for about 1 minute to allow the swab to completely absorb saliva.

Step 3: Put the swab on test machine.



## Miracle Technology Co.,Ltd.

Address: Room 301, Building A, No.18, Rongshuxia Industrial Zone, Tongxin Community, Baolong Street, Longgang District, Shenzhen, China, 518116

Tel: +86-0755- 89616773 Fax: +86-755-89616775

Web: [www.mantacc.com](http://www.mantacc.com)

EC REP

## Share Info Consultant Service LLC Repräsentanzbüro

Address: Heerdter Lohweg 83, 40549 Düsseldorf

Tel: 0049 1767 0057022

Fax: 0049 1767 0057022

Emai: [info@share-info.cn](mailto:info@share-info.cn)

DIMDI NO.: DE/0000047946

CE Logo

