

◆◆ Enteral Syringe

The Enteral syringe is used for delivering drug or food to oral or enteral

◆◆ Structure and composition:



- ◆ The Enteral syringe are assembled by Barrel, Plunger, Piston. All parts and material for this product meet medical requirements, after sterilized by ETO.

◆◆ Material:

- ◆ PP, Isoprene rubber, Silicone oil

◆◆ Specifications:

- ◆ 1ml,3ml,5ml,10ml,20ml,60ml,100ml



◆◆ Oral Syringe

Oral syringe for single use are intended for oral administration of drugs.

◆◆ Structure and composition:

- ◆ The oral syringe are assembled by Barrel, Plunger, Piston. All parts and material for this product meet medical requirements, after sterilized by ETO.

◆◆ Material:

- ◆ PP, Isoprene rubber, Silicone oil

◆◆ Specifications:

- ◆ 1ml,3ml,5ml,10ml,20ml

