



Cannulas



Comparison between Regular Wall, Thin Wall and Extra Thin Wall



Long bevel



Short bevel

GAUGE	Size		Outer diameter		Inner diameter		
	Nominal		Min	Max	Regular wall min	Thin wall min	Extra thin wall min
31G	0.25		0.254	0.267	0.114	0.125	
30G	0.30		0.298	0.320	0.133	0.165	
29G	0.33		0.324	0.351	0.133	0.190	
28G	0.36		0.349	0.370	0.133	0.190	
27G	0.40		0.400	0.420	0.184	0.241	
26G	0.45		0.440	0.470	0.232	0.292	
25G	0.50		0.500	0.530	0.232	0.292	
24G	0.55		0.550	0.580	0.280	0.343	
23G	0.60		0.600	0.650	0.317	0.370	0.460
22G	0.70		0.698	0.730	0.390	0.440	0.522
21G	0.80		0.800	0.830	0.490	0.547	0.610
20G	0.90		0.860	0.920	0.560	0.635	0.687
19G	1.10		1.030	1.100	0.648	0.750	0.850
18G	1.20		1.200	1.300	0.790	0.910	1.041
17G	1.40		1.400	1.510	0.950	1.156	1.244
16G	1.60		1.600	1.690	1.100	1.283	1.390
15G	1.80		1.750	1.900	1.300	1.460	1.560
14G	2.10		1.950	2.150	1.500	1.600	1.727
13G	2.40		2.300	2.500	1.700	1.956	
12G	2.70		2.650	2.850	0.950	2.235	
11G	3.00		2.950	3.150	2.200	2.464	
10G	3.40		3.300	3.500	2.500	2.819	

