-PRODUCTS-

http://www.mekonneedle.com

**** Enteral Syringe**

The Enteral syringe is used for delivering drug or food to oral or enteral

* Structure and composition:

 The Enteral syringe are assembled by Barrel, Plunger, Piston. All parts and material for this product meet medical requirements, after sterilized by ETO.


* Material:

• PP, Isoprene rubber, Silicone oil

* * Specifications:

1ml,3ml,5ml,10ml,20ml,60ml,100ml

MSK



** Oral Syringe

Oral syringe for single use are intended for oral adminstration of drugs.

* * Structure and composition:

 The oral syringe are assembled by Barrel, Plunger, Piston. All parts and material for this product meet medical requirements, after sterilized by ETO.

* Material:

PP, Isoprene rubber, Silicone oil

Specifications:

1ml,3ml,5ml,10ml,20ml

