



CHG Swabsticks

Products Descriptions

CHG Swabsticks are pre-saturated disinfectant swabs which have been applied for skin antiseptis preparation, they are made of medical grade 100PPI open-cell polyurethane foam & white ABS handle. It is with formulation of 2% chlorhexidine gluconate (CHG) and 70% isopropyl alcohol (IPA) to meet the strict clinical requirements of FDA, thus it would kill more bacteria than traditional Iodophor Swabs & Alcohol Swabs. CHG Swabsticks have been widely used for Hospitals, CDC (Center for Disease Control and Prevention), INS (Infusion Nurses Society), SHEA (Society for Healthcare Epidemiology of America) and IDSA (Infectious Diseases Society of America). Please be advised that customization of solution formulation & size of Swabstick is available.

Products Features

- Convenient for using and two-year shelf life.
- Latex-free and does not contain any DEHP.
- Persistent, friction scrub and aseptic technique.
- Active in protein-rich biomaterials and rapid-acting.
- Broad spectrum: MRSA, VRE, CRE, Clostridium difficile and Acinetobacter.

Price List of CHG Swabs

Miraclean Part No.	Product Descriptions	Prices Range (usd)			Packing Information
		5,000-10,000	10,000-50,000	50,000-200,000	
MIP-707F	MIP-707F Disinfection Swab with Rectangular Sponge Head	\$0.16	\$0.14	\$0.12	50pcs/box, 40boxes/case, the case dimension is 505 ×342 ×435 mm.
MIP-708F	MIP-708F Disinfection Swab with Circular Sponge Head	\$0.16	\$0.14	\$0.12	50pcs/box, 40boxes/case, the case dimension is 505 ×342 ×435 mm.

Products Specifications

Item No.	Total length	Head length	Head width	Head thickness	Handle length	Handle width	Handle thickness
MIP-707F	131 mm	25 mm	15 mm	8.8 mm	106 mm	6.6 mm	3.2 mm
Solution formulation, 2% Chlorhexidine gluconate and 70% Isopropyl alcohol							
Packing information, 100pcs /box, 30boxes / case, the case dimension is 480 ×320 ×250 mm							



Item No.	Total length	Head length	Head width	Head thickness	Handle length	Handle width	Handle thickness
MIP-708F	132 mm	26.5 mm	19 mm	8.5 mm	105.5 mm	6.6 mm	3.2 mm
Solution formulation, 2% Chlorhexidine gluconate and 70% Isopropyl alcohol							
Packing information, 100pcs / box, 30boxes / case, the case dimension is 380 ×340 ×320 mm							



Sterilization Process

We take Gama Ray sterilization process for CHG Swabsticks if customers have specific requirement.